HAYER & BOECKER

DIE DRAHTWEBER

HAYER USC 200 Multi

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Index

		<i>Page</i>
1.	<i>Important safety instructions</i>	12
2.	<i>Quantity delivered</i>	13
3.	<i>General operating instructions</i>	13
4.	<i>Starting the unit</i>	14
5.	<i>Further recommendations</i>	17
6.	<i>Technical Data</i>	18
7.	<i>Guarantee</i>	19
	<i>Annex: EC-Attestation of conformity</i>	

1. Important safety instructions

- After removing the packing, check the unit's casing and control elements for **any sign of damage**. If you should find any damage, do **not** connect the unit to the electrical power supply, but report the damage immediately and in writing to the carrier and to HAVER & BOECKER. Save the original packing.
- Place the unit in a **horizontal position** on a dry, solid base!
- Plug the unit into **grounded sockets** only!
- Replace defective parts **only** with original parts!
- Do not clean any living beings (plants or animals) in the unit!
- Ultrasonic cleaning produces cavitation noises. If working within 2 meters of the unit for longer periods of time, be sure to wear **ear plugs** to avoid damage to your hearing!
- To transport the **empty** unit, use the two fixed handles.

2. Quantity delivered

- 1 USC -unit
- 1 cover
- 1 power cable
- 1 run-off hose
- 1 sieve-holder USC 200 Multi
- 5 liter bottle of USC-cleaning fluid for USC 200 Multi

Observe the operating instruction when using HAVER Test Sieve Cleaners. If used improperly, the manufacturer will assume no responsibility for safety or proper working order of the unit.

For any after-sales services please contact the following address:

**HAVER & BOECKER, Wire Weaving Division, Ennigerloher Straße 64,
D-59302 OELDE, Phone: +49-2522-30330 EXPORT;
Fax: +49-2522-30404; Email: PA@haverboecker.com**

3. General operating instructions

- **Never** run the unit without liquid.
- Do **not** fill the stainless steel oscillation tank with: combustible liquids (e.g. benzene, solvents) or chemicals which contain or release halogen ions, such as some disinfectants, dish-washing detergents, household cleaners.
- Use aggressive cleaning agents such as acids or saline solutions **only** in inset baskets.
- Do **not** put anything on the bottom of the oscillation tank, use the sieve rack for the cleaning of test sieves.
- Do **not** put your hand into the oscillation liquid during the cleaning process.
- When running the unit for a long period of time, the ultrasonic energy in the oscillation tank heats up the cleaning liquid. Check the temperature when cleaning temperature-sensitive parts.

**HAYER-Test Sieve Cleaners are RFI-tested and conform to CE-mark
Limits as per EN 50081-1 and 50082-1**

4. Starting the unit

4.1 Cleaning Advice

During each cleaning process be sure that the test sieves are completely covered with cleaning fluid. The cleaning process normally takes place directly in the oscillation tank.


Indirect cleaning in inset baskets or in an insert tank is recommended for:

- the use of chemically aggressive cleaning liquids such as acids etc.
- the cleaning of chemically aggressive soilings
- the cleaning of abrasive soilings such as polish paste, resin or sand.

Principle of ultrasonic cleaning

PZT (lead zirconate titanate) oscillators underneath the tank convert electrical energy into mechanical vibration. HAVER Test Sieve Cleaners cause the cleaning liquid to oscillate at 35 kHz. In doing so, tiny vacuum bubbles form which then implode (cavitation) This cavitation principle removes dirt and grease from the deepest pores - "electronic scrubbing"- at all points that come into contact with the cleaning liquid.

Direct cleaning in the oscillation tank

- The cleaning liquid should fill 2/3 of the oscillation tank. If the filling height is lower, the unit may be damaged.
- The cleaning process should only take place with an inset basket or unit holder.
- Auxiliary parts should not be placed on the bottom of the oscillation tank.

When inserting the goods to be cleaned it is important to note that

- the most soiled side faces down
- the sieves are not stacked together, otherwise the ultrasonic action will be absorbed


- parts should not touch each other; they can be cleaned in holders or baskets preventing damage
- there should be no air bubbles in cavities (e.g. pocket holes)
- temperature has to be controlled for temperature-sensitive parts.

Degassing of the cleaning fluid increases the cleaning efficacy. Cleaning liquids contain dissolved gases (e.g. oxygen). Irradiate cleaning liquids which have just been poured into the oscillation tank or have been left in it for a longer period of time apx. 5 to 15 minutes before starting the cleaning process. During the degassing process, the cavitation noise changes and the loud degassing noises decrease toward the end of the degassing process, the unit appears to run more quietly. Lower noise levels indicate the end of the degassing process, not a reduction in ultrasonic energy.

4.2 Switching the ultrasonic energy on and off

TIME OPERATION

- Turn the knob to the right to set the appropriate time of cleaning → **the green pilot lamp lights up.** The timer switches off **automatically**.
- You can shorten the cleaning time or switch the unit off by turning the knob in the opposite direction.

CONTINUOUS OPERATION (∞ position)

- Turn the knob to the **left** until **it clicks into position**.
- The unit does not switch off automatically, turn the knob back to "0" in order to switch the unit off.


The unit may remain plugged in when it is not running. The power supply is interrupted when pulling out the plug from the mains supply.


4.3 Cleaning time

In principle, the cleaning time must be kept as short as possible. Depending on the soiling degree and the cleaning liquid (follow the instructions on the label of cleaning liquid) test sieves should be cleaned for approx. 1 to 15 minutes.

4.4 Emptying the oscillation tank

Do not put the unit into the sink. Turn the lever of the ball valve in the direction of outflow to open.

Mounting instructions for ball valve G ½

1. Wind the attached white PTFE band around the screw neck 2. Wind in clock-wise direction.
2. Wind with the PTFE band around thread of hose attachment 3.
3. Screw the hose attachment 3 into the outlet of the ball valve 4. Turn in a clock-wise direction using spanner SW 22.
4. Screw the ball valve onto the screw neck, turning in a clock-wise direction and tighten (about 5 turns). **Attention:** Do not over-tighten ball valve. Do not unscrew it as this will damage the seal created by the PTFE band.
5. Adjust the locking nut 1 with a SW 32 spanner to ensure that the ball valve is in an upright position.
6. **To close the ball valve** turn the handle 5 to the left.

5. Further Recommendations

5.1 Cleaning means

- USC-cleaning fluid has proved most reliable for a deep and thorough cleaning. Please follow the **EU-safety-sheet** that is added to the cleaning fluid.
Dilute the cleaning fluid in water 10 : 1 (1 part of cleaning fluid, 10 parts of water)
- A soiled bottom of the oscillation tank or insert baskets (e.g. abrasives, polish other sedimentations) reduces the performance. Therefore, rinse out the oscillation tank and clean it, replace the cleaning liquid.
- Cleaning- or contact fluid in the ultrasonic oscillation tank should not exceed a maximum temperature of 100°C.
- Distilled or deionized water is not suitable for an ultrasonic cleaning directly in the oscillation tank. It has no cavitation germs to produce an even ultrasonic field. An irregular cavitation causes selective material wear of the products to be cleaned and the oscillation tank.
- Only use the USC-cleaning fluid in the oscillation tank.
- Distilled or deionized water may be used in inset baskets as cleaning fluid or to dilute other cleaning fluids.
- Take care when using aggressive cleaning agents in inset baskets, avoid splashes into the contact liquid or on the stainless steel surface. If this happens, exchange the contact fluid immediately, clean and wipe the surfaces dry.

The guarantee for damaged units or instruments is void if improper cleaning- or disinfection chemicals have been used.

5.2 Cleaning and care of the unit

- Be sure to unplug the unit before cleaning.
- Do not rinse or dip the unit in water, but humidly wipe its in- and outside.
- Rinse out the stainless steel tank frequently and wipe dry.
- If rings should appear on the stainless steel surface of the oscillation tank after a longer period of operation, remove them with a commercially available stainless steel cleaner without scouring additives. Do not use steel wool, scrapers or graters.
- Metal parts remaining on the stainless steel surface, as well as particles of rust, can penetrate the passive protective layer of the stainless steel, thereby "activating" it and causing it to rust. This rust can then pit the stainless steel surface. Therefore, remove all remaining metal parts; remove minor rust spots immediately with a soft cloth and commercially available stainless steel cleaner.

Problem Solving

- **The unit's oscillation is not as strong, uniform or quiet as it should be:**
 - Has the liquid been properly degassed? → **irradiate for 15 minutes!**
 - Irregular noises (wobble) are not a defect → **change the level of the liquid!**
 - Damages on the bottom of the tank caused by erosion.
 - Overloaded with products to be cleaned? Remove some of them.
- Unit is equipped with an overvoltage suppressor.
- **The unit only should be repaired by authorized experts!**

6. Technical Data

Operation voltage	220-240V; 50-60Hz, N, PE (available with special voltage 110V)
HF-Frequency	35kHz
Oscillation tank (inside)	500 x 300 x 300mm
Charge	40 l
Outlet	Ball valve G ¹ / ₂
Outside dimensions	560 x 360 x 500mm
Weight	27 kg
HF-constant maximum output	2 x 1000W/Per.
Power consumption	4,5 A
Fuses	F6A

Environmental conditions according to EN 61 010-1 (IEC 1010-1)

Degree of soiling	2 according to IEC 664
Overvoltage category	II
allowable surrounding temperature	5° to 40°
allowable relative humidity up to 31 °C	80%
allowable relative humidity up to 40 °C	50%
defrosting not allowable	

7. Guarantee

*All functions of the sieving machine are checked and adjusted before its delivery. The unit guarantee is for **24 months**.*

We are only liable for defects resulting from flaws or faults in material. The guarantee does not apply to damages caused by transport. Repairs resulting by improper handling or wear due to operation are excluded from guarantee. The guarantee is void, if interventions such as variations or improper repairs are effected which are outside our control.

Technical changes reserved

*HAYER & BOECKER
PARTICLE SIZE ANALYSIS*

*Ennigerloher Straße 64
59302 OELDE, Germany
Phone: 02522-30363 (GERMANY)
+49-2522-30330 (ABROAD)
Fax: +49-2522-30404
E-Mail: PA@haverboecker.com*

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IMDG: International Maritime Code of Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals`

Trade name	USC Cleaning Fluid	Revision:	17.09.2012
------------	--------------------	-----------	------------

1. Identification of the Substance / Mixture and of the Company / Undertaking

- **Product identifier**
- **Trade name: USC Cleaning Fluid**
- **Application of the substance / the preparation** Cleaning product
- **Manufacturer / Supplier:**
Haver & Boecker OHG
Ennigerloher Str. 64
D-59302 Oelde
Tel.: +49 2522 30-0
Fax: +49 2522 30-404
Mail: dw@haverboecker.com
- **Further information obtainable from:** Conservation department
- **Emergency telephone number:**
Informationszentrale gegen Vergiftungen
der Universität Bonn (GIZ-Bonn)
Adenauerallee 119
D-53113 Bonn

Tel: +49 228-19240
gizbn@mailen.meb.uni-bonn.de
www.meb.uni-bonn.de/giftzentrale/

2. Hazards Identification

- **Hazard description:**
- **Classification of the substance or mixture**
- **Classification according to Regulation (EC) No.1272/2008**

GHS05 Corrosion

Skin Corr. 1 B H314 Causes severe skin burns and eye damage.

- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**

Xi Irritant

R 36/38: Irritating to eyes and skin.

- **Information concerning particular hazards for human and environment:**
Preparation number 1999/45/EC in connection with substance directive 67/548/EEC.
The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.
- **Classification system:**
The classification is according to the latest editions of the EU-lists, and extended by company and literature data.
- **Label elements**
- **Labelling according to EU guidelines:**
Preparation number 1999/45/EC in connection with substance directive 67/548/EEC.
The product has been classified and marked in accordance with EU Directives / Ordinance on Hazardous Materials.
- **Code letter and hazard information of the product:**
Xi Irritant
- **Risk phrases:**
36/38 Causes eye and skin irritation.

Trade name	USC Cleaning Fluid	Revision:	17.09.2012
------------	--------------------	-----------	------------

- **Safety phrases:**
 - 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
 - 28 After contact with skin, wash immediately with plenty of water.
 - 36/37/39 Wear suitable protective clothing, gloves and eye / face protection.
 - 45 In case of accident or if you feel unwell, seek medical advice immediately.

3. Composition / Information on Ingredients

- **Chemical characterization:** Mixture
- **Description:** Mixture of substances listed below with nonhazardous additions.

Dangerous components:		
CAS: 68584-22-5	Benzenesulfonic acid, C10-16-alkyl derivs. C R34, Xn; R 22	2.5-10%
	Skin Corr. 1B ; H314 ; Acute Tox. 4 ; H302	
CAS: 69011-36-5 Polymer	Isotridecanol, ethoxyliert (6-9 EO) Xn R22 Xi R41	2.5-10%
	Eye Dam. 1, H318 ; Acute Tox. 4 ; H302	
CAS: 68585-34-2	Sulfohydroxypoly (oxy-1,2-ethandiy)C10-C16-alkyether,Na-Salze Xi; R 36/38	2.5-10%
	Skin Irrit. 2,H315; Eye Irrit .2, H319	
CAS: 34590-94-8 EINECS: 252-104-2	(2-methoxymethylethoxy)propanol Substance with a community workplace exposure limit.	2.5-10%
CAS: 77-92-9 EINECS:201-069-1	Citric acid Xi; R 36	2.5-10%
	Eye Irrit. 2 H319	

4. First Aid Measures

- **After inhalation:** In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:**
Immediately wash with water and soap and rinse thoroughly.
Immediately rinse with water.
- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** If symptoms persist consult doctor.
- **Information for the doctor:**
- **Most important symptoms and effects, both acute and delayed.**
No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed.**
No further relevant information available.

5. Fire-fighting Measures

- **Extinguishing agents:**
- **Suitable extinguishing agents:**
CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Protective equipment:** No special measures required.

7. Handling and Storage

- **Handling:**
- **Precautions for safe-handling:** Keep receptacles tightly sealed.
- **Information about fire- and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.

Trade name	USC Cleaning Fluid	Revision:	17.09.2012
------------	--------------------	-----------	------------

8. Exposure Controls / personal Protection

- **Additional information about design of technical facilities:** No further data; see item 7.

Ingredients with limit values that require monitoring at the workplace:

34590-94-8 (2-methoxymethylethoxy)propanol

WEL	Long-term value: 308 mg/m ³ , 50 ml/m ³
	Sk

- **Additional information:** The lists valid during the making were used as basis.

- **Personal protective equipment:**
- **General protective and hygienic measures:**
Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.
- **Respiratory protection:** Not required.
- **Protection of hands:**

Protective gloves

The glove material has to be impermeable and resistant to the product / the substance / the preparation. Due to missing tests no recommendation to the glove material can be given for the product/the preparation/ the chemical mixture.

- **Material of gloves:**
- **Protection of hands: Protective gloves made from plastics or rubber (EN 374).**
Material: Butyl rubber

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material cannot be calculated in advance and has therefore to be checked prior to the application.

- **Material thickness:** -0,5 mm
- **Penetration time of glove material:** >8 h
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **Eye protection:**

Tightly sealed goggles

9. Physical and chemical Properties

General information:

Form:	fluid
Colour:	yellowish - colourless
Odour:	characteristic
Odour threshold:	undetermined

Change in condition:

Melting point/Melting range:	undetermined
Boiling point/Boiling range:	100 °C
Flash point:	not applicable
Ignition temperature:	270°C
Decomposition temperature:	undetermined
Self-igniting:	Product is not self-igniting.

(Contd. on page -4-)

Trade name	USC Cleaning Fluid	Revision:	17.09.2012
------------	--------------------	-----------	------------

Danger of explosion:	Product does not present an explosion hazard.
Vapour pressure at 20°C	23 hPa
Density at 20°C :	ca. 1.01 g/cm ³
Solubility in/Miscibility with water:	Fully miscible.
pH-value at 20°C :	ca. < 4

10. Stability and Reactivity

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11. Toxicological Information

- **Acute toxicity:**
- **Primary irritant effect:**
- **On the skin:** Irritant to skin and mucous membranes.
- **On the eye:** Irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**
The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:
Irritant

12. Ecological Information

General notes:

Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water
Do not allow product to reach ground water, water course or sewage system.
Do not dump (unfiltered) product as concentrate into sewage water.

13. Disposal Considerations

- **Product:**
- **Recommendation:**
Must not be disposed together with household garbage. Do not allow product to reach sewage system.
- **European Waste Catalogue:**
According to EWC regulation, the allocation of waste code numbers is to be carried out specific to the industry and the process. No waste code according to the European Waste Catalogue (EWC) can be assigned for this product, as the intended use dictates the assignment. The waste code is established in consultation with the regional waste disposer.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

14. Transport Information

• UN-Number	
• ADR/ADN IMDG, IATA:	not applicable
• UN proper shipping name	
• ADR,ADN,IMDG,IATA:	not applicable
• Transport hazard class(es)	
• ADR, ADN, IMDG, IATA	
• Class	not applicable
• Environmental hazards:	
• Marine pollutant:	no
• Special precautions for user:	not applicable
• Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	not applicable

Trade name	USC Cleaning Fluid	Revision:	17.09.2012
------------	--------------------	-----------	------------

15. Regulatory Information

- Safety, health and environmental regulations/legislation specific for the substance or mixture.
- National regulations:
- Technical instructions (air):

Class	Share in %
NK	2.5-10

- Water hazard class: Water hazard class 2 (self-assessment): hazardous for water.
- Chemical safety assessment: A chemical safety assessment has not been carried out.

16. Other Information

These data are based on our current level of knowledge and experience. However, they shall not constitute a guarantee for any specific product features and justify no contractual legal relationship.
Department issuing data sheet: Conservation department

Abbreviations and acronyms:

ADR: Accord europeen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code of Dangerous Goods
IATA: International Air Transport Association
GHS: Globally Harmonized System of Classification and Labelling of Chemicals`